

# EARLSTON'S DEFIBRILLATORS

## WHERE ARE THEY?

1. EARLSTON PRIMARY SCHOOL  
On the wall next to the main school entrance
2. RED LION HOTEL, EARLSTON  
On the gable wall on facing the pharmacy
3. EARLSTON RUGBY CLUB  
On the gable wall nearest the road

## WHO CAN USE IT?

In an emergency anyone can, you don't need to be trained.

## HOW DO I ACCESS IT?

The box has no lock or security device, just open the door

## HOW DOES IT WORK?

Just follow the instructions and the defibrillator will tell you what to do.

**CALL EMERGENCY  
SERVICES ON 999**



# PERFORMING A RESCUE

1




Is patient responsive?

2



Call emergency services!

3



Get AED.

4



Place AED on ground.  
Open lid.

5



Remove clothing.  
Expose chest.

6



Attach pads.

7



Do not touch the patient.

8



If required, AED will deliver  
shock automatically.

9



If present, place CPR device.

10



Follow AED prompts.  
Give compressions.

11



If instructed, give breaths.

## Automatic AED



Cardiac Science and the Shielded Heart logo are trademarks of Cardiac Science Corporation.

Copyright © 2012 All Rights Reserved.

N7 W22025 Johnson Drive, Waukesha, WI 53186 USA  
[www.cardiacscience.com](http://www.cardiacscience.com)



MDSS GmbH  
Schiffgraben 41  
D-30175 Hannover  
Germany



70-01092-01r E 0086

